REPORT	DOCUM	ENTATION PAGE			Form Approved OMB No. 0704-0188
The public reporting burden for this collection of gathering and maintaining the data needed, and conformation, including suggestions for reducing that notwithstanding any other provision of law, control number.  PLEASE DO NOT RETURN YOUR FO	ompleting and ne burden, to no person si	d reviewing the collection of infor the Department of Defense, Ex- nall be subject to any penalty for	mation, Send com ecutive Services a r failing to comply	luding the tir ments regard nd Communic with a collec	me for reviewing instructions, searching existing data sources, ding this burden estimate or any other aspect of this collection of cations Directorate (0704-0188). Respondents should be aware ction of information if it does not display a currently valid OMB
1. REPORT DATE (DD-MM-YYYY) 13-03-2003		ORT TYPE Final Technical			3. DATES COVERED (From - To) 1 Jan 2002 31 Dec 2003
4. TITLE AND SUBTITLE	<u> </u>	· · · · · · · · · · · · · · · · · · ·		5a. CO	NTRACT NUMBER
The Measurement of Trace Metals	in JP-8 b	y ICP-OES, Presentati	on 1990-8	5b. GR/	ANT NUMBER
				5c. PRO	OGRAM ELEMENT NUMBER 63112F
6. AUTHOR(S)				5d. PRC	DJECT NUMBER
Timothy J. Shelley, Ph.D., Air For		rch Laboratory, Airbase	e		
Technologies Division, Tyndall Al	PB FL			5e. TAS	SK NUMBER
Christian Voelkl, Ph.D., Bundeswe	hr Resea	rch Institute for Materi	als,		
Explosives, Fuels and Lubricants,	Erding, C	Germany		5f. WO	RK UNIT NUMBER
					4918L41C
7. PERFORMING ORGANIZATION NA	ME(S) A	ND ADDRESS(ES)			8. PERFORMING ORGANIZATION
Det 2, Air Force Research Laborat	ory				REPORT NUMBER
139 Barnes Drive, Suite 2					AFRL-ML-TY-TP-2004-4536
Tyndall AFB FL 32403					
9. SPONSORING/MONITORING AGE	NCY NAM	E(S) AND ADDRESS(ES)			10. SPONSOR/MONITOR'S ACRONYM(S)
Det 2, Air Force Research Laborat					
139 Barnes Drive, Suite 2					
Tyndall AFB FL 32403					11. SPONSOR/MONITOR'S REPORT
12. DISTRIBUTION/AVAILABILITY ST	ATEMEN	Г	-		AAA/AAAA AAA
Distribution Unlimited; Approved	for Publi	c Release		*: 	20040802 081
13. SUPPLEMENTARY NOTES				<del></del>	
	2003 (Pi	ttsburgh Conference or	Analytical (	Chemistry	y and Applied Spectroscopy) in Orlando, FL
14. ABSTRACT	<b>4</b>		. T.,	TD 0	
JP-8 is the preferred fuel for Air Fe	orce trans	sport and fighter aircrate, boiling range, and the	t. In general ermal stabilit	, JP-8 acc v. Verv l	ceptance specifications focus on physical little work has been performed to identify
and measure the metals that may b	e present	in JP-8 at an operating	airbase. The	ere is an i	interest in the catalytic processing of JP-8
					s associated with airbase operations. The
goal of this project was to develop discuss some of the challenges und	an analy	tical method for the det	termination o	f metals l	by direct aspiration ICP/OES. We will
discuss some or the chancinges und	overed d	aring method developing	nem and pro-	ride resur	is generated by the method.
15. SUBJECT TERMS					
fuels, JP-8, catalyst, inorganic che	mistry				
16. SECURITY CLASSIFICATION OF:		17. LIMITATION OF	18. NUMBER	19a. NA	ME OF RESPONSIBLE PERSON
a. REPORT   b. ABSTRACT   c. Th		ABSTRACT	OF		y J. Shelley, Ph.D.
		7777	PAGES	10h TEL	EDUONE MUMPED //palvide area godel

U

U

U

UU

16

19b. TELEPHONE NUMBER (Include area code) 850-283-3605



# The Measurement of Trace Metals in JP-8 by ICP-OES

Fimothy J. Shelley, Ph.D. Air Force Research Laboratory, Airbase Technologies Division, Tyndall AFB, FL Christian Voelkl, Ph.D., Bundeswehr Research Institute for Materials, Explosives, Fuels and Lubricants, Erding Germany









## ≯Organizational Background

- ✓ Description
- ➤ Resources Profile
- ➤ Research Thrust Areas


# > Research Areas and Programs

- ≻Robotics Research
- >Englineening Mechanics
  - Explosives Effects
- 4 Expalorative Ops & Stuppoulit

March 13, 2003



### Air Force Research Laborator TECHNOLOGY DIRECTORALES



March 13, 2003



### Air Force Research Laboratory LOCATIONS



### Space Vehicles Sensors AFOSR (Wash DC) Headquarters Air Vehicles Materials & Mfg Propulsion Sensors

Information

ROME

Sensors

**EDWARDS** 

Propulsion

March 13, 2003

Effectiveness

Human

MESA

Approved for public release

Materials & Mfg

**LYNDAL** 

Munitions

EGLIN

**Effectiveness** 

Human

BROOKS

**Directed Energy** 

Space Vehicles

KIRTLAND

4

Human Effectiveness

Information



# 

Aliribasse Technicologijes Divisioni MLO

Operations Support Branch (MLQO)

Weapons Systems Logistics Branch (MLQL)

> Deployed Base Systems Branch (MLQD)

Approved for public release

Force Protection Branch (MLQF)

March 13, 2003




# AFRL - Tyndall AFB

Force Protection Branch MLQF Airbase Technologies Division



# > Research Areas and Programs

- ➤ Robotics Research
- ➤ Engineering Mechanics
- ➤ Explosives Effects
- YEXPlosive Ops & Support
- ➤ Chreimical/Biologologi Delense




March 13, 2003


Approved for public release



# Jeployable Fuel Cell Reforme







■ 50% Reduction in Power Deployment Airlift

CLEAN

HEAT AND WATER

- Increase in MTBF from 500 to 2200 hrs
  - Savings of 1800 gallons of fuel/day (1100 man deployment)
    - Reduced Acoustic / Thermal / Environmental Emissions

POWER

POWER


PROCESSOR


AIR




#### **ICP System**







- Thermo Iris Advantage
  Dual View (Axial)
  Simultaneous ICP
  175-800 nm
- GE ABC Fully Demountable Torch
- GE MicroMist Nebulizer







#### **Objective**

- ➤ Minimal Sample Prep
  - > Sensitive
- ➤Metals, Phosphorus, Sulfur

### ➤ Challenges➤ Prior Work➤ Blanks➤ Matrix Effects

- >Standards
- >Sampling/Storage/Stabilin

March 13, 2003





#### ➤ Blank Solutions

JP-8: S; Cu, Si, P; some Al, Fe


Kerosenes: Al, P

**Conostan Premisolv** 

#### Solvent Uptake Rate

Neb flow	pump	Premisolv	Kerosene	JP-8
Ipm	rpm	mg/min	mg/min	mg/min
0.55	80	41.9	45.4	48.7
0.5		35.1	39.6	43.5
0.55	70	40.0	46.0	48.6
0.5		35.7	40.8	45.2
0.55	09	36.5	39.3	44.6
0.5		30.9	36.3	42.6
0.55	50	29.3	34.4	41.0
0.5	50	28.5	32.6	37.7







### Analytical Conditions



Solveng 80:20 mixture (wolve) Prendsolv Xylene

Nebulizer: GE MicroMist, 400 uL/min

Injector tube: 1mm I.D.

Tubing: Glass Expansion, Viton, orange/yellow

Pump seed: 60 rpm


Nebulizer Flow: 0.53 I/min (.50 - .55)


Sample Uptake: 40 mg/min


Internal Std: Sc2554, Sc3353; 800 ppb



#### Results







Sample Storage:

Glass: Sorption; Si

Polyethylene: Phosphorus Leaching

Working Standards: daily prep





### <u>Acknowledgements</u>

Christian Voelkl, Ph.D., Bundeswehr Research Institute for Materials, Explosives, Fuels and Lubricants, Erding Germany

Jerry Dulude, Glass Expansion Inc.

Christia Hodges, Thomas Dicorda Corp.